

# BabySensor QUICK GUIDE











# CONNECT THE BASE AND BRACELET

- Place the base and sensor on charge
- Remember to remove the plastic under the sensor
- Fully charge them the first time



- Make sure your phone is connected to a network that supports 2.4 GhZ
- Most networks support both 2.4 GhZ and 5 GhZ, or are split into two different networks
- Download the app from Apple Store / Play Store
- Follow the instructions in the app
- Test that you receive the first heart rate measurement

#### **CONNECT THE CAMERA**

- Plug in the camera
- The camera announces it is ready for setup
- Check that your phone is on the same network as the base
- Follow the instructions in the app, under Settings / Live stream



 Choose settings for sound, light mode, etc.




#### **FAMILY SHARING**

- Check that everything works for you before sharing
- Go to Settings / Share with family member in the app. Follow the instructions



- The person you share with receives an invitation, downloads the app, and logs in with the account you created for them
- Each user can choose their own notification settings







#### **CONNECT THE BASE AND BRACELET**



Cannot connect the base to WiFi. Message says wrong network or password.



Make sure the password is correct and that the network supports 2.4 GhZ.

# Solution

- 2.4 GhZ networks have better range but are slower than 5 GhZ
- 5 GhZ networks often have "5" in the name
- Some networks can be difficult to connect to. Ensure good coverage.
- If you have your own router, set it to send 2.4 GhZ. Remember to refresh the router after changes.



2



App does not find the camera for setup or use



Make sure the camera is connected to power (also for reset).

# Solution

- Tilt the black part of the camera upwards
- At the bottom, you'll see a port and a small black button to the left
- Hold the button until the camera says "Factory reset"
- Follow the app instructions under Settings / Live stream

#### **FAMILY SHARING**

Thrown out of the app or not receiving alerts as expected

# Check

Check "Member access" to ensure sharing is enabled and used.

# Solution

- Go to Settings / Member access
- Delete the invited user
- Ask the invited user to delete the app
- Make sure everything works for you
- Share again via Settings / Share with family member







#### **REFRESH**



# Problem

The app is frozen or functions are not working



### ← Sjekk

Make sure you have the latest version of the app and desired functions are enabled.



### Løsning

- Close and restart the app
- Restart the base station (use the power button on the front or unplug/replug power cable)

#### **EASY RESET**



# Problem

Several functions missing or not working



# Sjekk

Make sure you have the latest version of the app and desired functions are enabled.



#### Løsning

- Go to Settings / your name (top of next page) / Log out
- Delete the app
- Download the new app
- Log in with your existing user
- Restart the base station

#### **FULL RESET**



# Problem

No contact between app and base station, not receiving bracelet data



## Sjekk

Make sure the base is connected to WiFi or 4G. Check that the bracelet blinks green or red when against skin.

# ( Løsning

- Go to Settings / your name (top of next page) / Deregister account
- Delete the app
- Download the app again
- Set up the system with a new user
- Make sure everything works
- Share again with family member